

## Excel Spreadsheet Program for Viewing echoPlaque Measurements

One of the key reporting tools used in echoPlaque is the powerful DAT Summary Program. After performing analysis and generating measurements in echoPlaque, your results can be exported and displayed using the DAT Summary Program. Included with each copy of echoPlaque, the DAT Summary Program is an extensive, interactive Microsoft Excel based spreadsheet application that contains frame by frame measurements, region measurements, summaries, averages, and more.

Complete measurements organized over regions and frames:

Version: 3 single treatment range		Treatment and Two Reference Regions							
Summary from C:\echovis\data\Remodeling\IndexCompact.dat									
Frame	Pos	MaxLD	MinLD	AvgLD	LArea	MaxVD	MinVD	AvgVD	
192	6.40								
221	7.37								
	1.00	2.99	2.29	2.65	5.25	5.97	4.95	5.40	
222	7.49								
393	13.10								
	5.73	3.26	2.52	2.65	6.46	5.97	4.75	5.35	
394	13.13								
423	14.10								
	1.00	3.64	2.24	2.89	6.65	5.46	4.22	4.75	
Frame	Pos	MaxLD	MinLD	AvgLD	LArea	MaxVD	MinVD	AvgVD	
192	6.40	3.53	2.51	3.01	1.00	7.51	5.80	6.40	
193	6.43	3.45	2.50	2.98	6.90	7.16	5.78	6.31	
194	6.47	3.45	2.47	2.95	6.74	7.06	5.72	6.24	
195	6.50	3.43	2.47	2.94	6.71	6.88	5.69	6.17	
196	6.53	3.32	2.44	2.90	6.54	6.84	5.64	6.11	
197	6.57	3.27	2.42	2.89	6.45	6.73	5.61	6.03	
198	6.60	3.26	2.42	2.86	6.34	6.62	5.53	5.95	

Summation of key measurement results:

