


## Review and Analysis of iMap™ data using INDEC Medical software

INDEC's Enhancer software is now fully compatible with your Boston Scientific iMap™ tissue classification datasets. Review and analysis of iMap data can now be performed with INDEC's user-friendly software that can be added to an existing INDEC system or purchased as standalone software for your desktop or laptop computer\*.

- Load grayscale pullbacks with full color iMap data and contours
- Multiple load modes (iMap frames only or full pullback with iMap)
- Create new or edit existing contours
- Choose which tissue types are included in classification
- Confidence Level threshold to refine results
- Export measurement results to reports
- Save images and animations for presentations
- Compatible with INDEC's Liberator storage software - query, retrieve, and store your iMap data


### Display Boston Scientific iMap images, contours, and measurement results:






INDEC Systems, Inc. 2210 Martin Ave. Santa Clara, CA 95050 Tel: 1-408-986-1600

[www.indecsystems.com](http://www.indecsystems.com)


Edit or create new measurements. iMap coloring and measurements are automatically updated:




Additional features to improve and enhance analysis:



Export results to Excel Spreadsheets or create images and animations for presentations:



Source Frame Position	Lumen Area	Vessel Area	Fibrotic %	Fibrotic Conf.	Lipidic %	Lipidic Conf.
25	0	20.39	29.01	74%	81%	4%
56	0.91	17.41	33.28	72%	79%	5%
86	1.4	15.32	23.3	83%	85%	4%
116	1.9	13.93	22.05	84%	85%	4%
146	2.39	14.56	21.67	78%	82%	4%
176	2.89	13.75	22.06	77%	83%	4%
206	3.38	15.14	22.74	69%	81%	3%
235	3.86	14.23	16.87	84%	88%	1%
265	4.36	14.86	23.2	65%	80%	4%
295	4.85	14.91	22.81	74%	79%	4%
324	5.33	15.67	22.47	76%	79%	7%
354	5.82	14.4	22.24	83%	82%	4%
384	6.32	14.52	24.2	82%	83%	5%
414	6.81	15.94	22.49	89%	83%	2%
443	7.29	13.66	19.87	91%	87%	5%
473	7.79	12.54	19.71	88%	83%	2%
502	8.27	14.41	20.9	78%	80%	3%



INDEC Systems, Inc. 2210 Martin Ave. Santa Clara, CA 95050 Tel: 1-408-986-1600

[www.indecsystems.com](http://www.indecsystems.com)

iMap is a trademark of Boston Scientific Corporation and is not available in the US  
 \*Computers running Enhancer software must meet minimum system requirements